

THE FIBRICHECK APP



FibriCheck is a smartphone application that is capable of recording your heart rate and heart rhythm. It does this by making a connection between your finger and the smartphone camera. In the instructions below you'll find a short explanation how-to-get-started and how-to-use the application.


Step 1.

- Go to the “play store” or the “Apple App Store” on your smartphone.
- Search the app ‘FibriCheck’ and download the application.



Step 2.

After the first-login you'll need to create an account, fill in your own personal data and the email address you provided.



Step 3.

- To make a registration of the heart rhythm, click “start”
- place your index finger “gently” on the camera of your smartphone (you do not need to cover the flashlight, this can become hot!). The registrations will start automatically and will take 60 seconds. Afterwards the data is automatically sent to your caregiver.

